

CLIENT DETAILS

Contact Luke Dewson  
 Client NAMOI COTTON LIMITED  
 Address PO BOX 1333  
 TOOWOOMBA QLD 4350

Telephone 61 7 46316100  
 Facsimile 61 7 46316184  
 Email ldewson@namoicotton.com.au

Project **Merah North Cotton Gin Stormwater Discha**  
 Order Number **186485**  
 Samples 1

LABORATORY DETAILS

Manager Huong Crawford  
 Laboratory SGS Alexandria Environmental  
 Address Unit 16, 33 Maddox St  
 Alexandria NSW 2015


Telephone +61 2 8594 0400  
 Facsimile +61 2 8594 0499  
 Email au.environmental.sydney@sgs.com

SGS Reference **SE250482 R0**  
 Date Received 10 Jul 2023  
 Date Reported 17 Jul 2023

COMMENTS

Accredited for compliance with ISO/IEC 17025 - Testing. NATA accredited laboratory 2562(4354).

SIGNATORIES



Dong LIANG  
 Metals/Inorganics Team Leader

Sample Number SE250482.001  
 Sample Matrix Water  
 Sample Date 06 Jul 2023  
 Sample Name Point 18

Parameter	Units	LOR
-----------	-------	-----

**Anions by Ion Chromatography in Water Method: AN245 Tested: 13/7/2023**

Nitrate Nitrogen, NO <sub>3</sub> -N	mg/L	0.005	<b>0.56</b>
--------------------------------------	------	-------	-------------

**Nitrite in Water Method: AN277 Tested: 11/7/2023**

Nitrite Nitrogen, NO <sub>2</sub> as N	mg/L	0.005	<b>0.044</b>
Total Oxidised Nitrogen, NO <sub>x</sub> -N	mg/L	0.005	<b>0.61</b>

**TKN Kjeldahl Digestion by Discrete Analyser Method: AN292 Tested: 11/7/2023**

Total Kjeldahl Nitrogen	mg/L	0.05	<b>1.4</b>
Total Nitrogen (calc)	mg/L	0.05	<b>2.0</b>

**Total Phosphorus by Kjeldahl Digestion DA in Water Method: AN279/AN293(Sydney only) Tested: 11/7/2023**

Total Phosphorus (Kjeldahl Digestion) as P	mg/L	0.02	<b>0.49</b>
--	------	------	-------------

**pH in water Method: AN101 Tested: 10/7/2023**

pH**	No unit	-	<b>6.9</b>
------	---------	---	------------

**Conductivity and TDS by Calculation - Water Method: AN106 Tested: 10/7/2023**

Conductivity @ 25 C	µS/cm	2	<b>180</b>
---------------------	-------	---	------------

Sample Number	SE250482.001
Sample Matrix	Water
Sample Date	06 Jul 2023
Sample Name	Point 18

Parameter	Units	LOR
<b>Total and Volatile Suspended Solids (TSS / VSS) Method: AN114 Tested: 13/7/2023</b>		
Total Suspended Solids Dried at 103-105°C	mg/L	5
		<b>150</b>

<b>Oil and Grease in Water Method: AN185 Tested: 14/7/2023</b>		
Oil and Grease	mg/L	5
		<5

MB blank results are compared to the Limit of Reporting

LCS and MS spike recoveries are measured as the percentage of analyte recovered from the sample compared the the amount of analyte spiked into the sample.

DUP and MSD relative percent differences are measured against their original counterpart samples according to the formula : *the absolute difference of the two results divided by the average of the two results as a percentage*. Where the DUP RPD is 'NA' , the results are less than the LOR and thus the RPD is not applicable.

**Anions by Ion Chromatography in Water Method: ME-(AU)-[ENV]AN245**

Parameter	QC Reference	Units	LOR	MB	DUP %RPD	LCS %Recovery	MS %Recovery
Nitrate Nitrogen, NO3-N	LB285158	mg/L	0.005	<0.005	0 - 1%	98%	90%

**Conductivity and TDS by Calculation - Water Method: ME-(AU)-[ENV]AN106**

Parameter	QC Reference	Units	LOR	MB	DUP %RPD	LCS %Recovery
Conductivity @ 25 C	LB284918	µS/cm	2	<2	0%	94%

**Nitrite in Water Method: ME-(AU)-[ENV]AN277**

Parameter	QC Reference	Units	LOR	MB	DUP %RPD	LCS %Recovery	MS %Recovery
Nitrite Nitrogen, NO2 as N	LB284948	mg/L	0.005	<0.005	7%	96%	80%

**Oil and Grease in Water Method: ME-(AU)-[ENV]AN185**

Parameter	QC Reference	Units	LOR	MB	LCS %Recovery
Oil and Grease	LB285321	mg/L	5	<5	105%

**pH in water Method: ME-(AU)-[ENV]AN101**

Parameter	QC Reference	Units	LOR	LCS %Recovery
pH**	LB284918	No unit	-	100%

**TKN Kjeldahl Digestion by Discrete Analyser Method: ME-(AU)-[ENV]AN292**

Parameter	QC Reference	Units	LOR	DUP %RPD	MS %Recovery
Total Kjeldahl Nitrogen	LB284947	mg/L	0.05	5 - 15%	116%

MB blank results are compared to the Limit of Reporting

LCS and MS spike recoveries are measured as the percentage of analyte recovered from the sample compared the the amount of analyte spiked into the sample.

DUP and MSD relative percent differences are measured against their original counterpart samples according to the formula : *the absolute difference of the two results divided by the average of the two results as a percentage*. Where the DUP RPD is 'NA' , the results are less than the LOR and thus the RPD is not applicable.

**Total and Volatile Suspended Solids (TSS / VSS) Method: ME-(AU)-[ENV]AN114**

Parameter	QC Reference	Units	LOR	MB	DUP %RPD	LCS %Recovery
Total Suspended Solids Dried at 103-105°C	LB285194	mg/L	5	<5	1%	99%

**Total Phosphorus by Kjeldahl Digestion DA in Water Method: ME-(AU)-[ENV]AN279/AN293(Sydney only)**

Parameter	QC Reference	Units	LOR	MB	DUP %RPD	LCS %Recovery	MS %Recovery
Total Phosphorus (Kjeldahl Digestion) as P	LB284947	mg/L	0.02	<0.02	10 - 78%	106%	111%

METHOD

METHODOLOGY SUMMARY

AN101	pH in Soil Sludge Sediment and Water: pH is measured electrometrically using a combination electrode (glass plus reference electrode) and is calibrated against 3 buffers purchased commercially. For soils, an extract with water is made at a ratio of 1:5 and the pH determined and reported on the extract. Reference APHA 4500-H+.
AN106	Conductivity and TDS by Calculation: Conductivity is measured by meter with temperature compensation and is calibrated against a standard solution of potassium chloride. Conductivity is generally reported as $\mu\text{mhos/cm}$ or $\mu\text{S/cm}$ @ 25°C. For soils, an extract with water is made at a ratio of 1:5 and the EC determined and reported on the extract, or calculated back to the as-received sample. Total Dissolved Salts can be estimated from conductivity using a conversion factor, which for natural waters, is in the range 0.55 to 0.75. SGS use 0.6. Reference APHA 2510 B.
AN106	Salinity may be calculated in terms of NaCl from the sample conductivity. This assumes all soluble salts present, measured by the conductivity, are present as NaCl.
AN114	Total Suspended and Volatile Suspended Solids: The sample is homogenised by shaking and a known volume is filtered through a pre-weighed GF/C filter paper and washed well with deionised water. The filter paper is dried and reweighed. The TSS is the residue retained by the filter per unit volume of sample. Reference APHA 2540 D. Internal Reference AN114
AN185	Gravimetric Oil & Grease and Hydrocarbons: A known volume of sample is extracted using an organic solvent and the solvent layer with dissolved oils and greases is transferred to a pre-weighed beaker. The solvent is evaporated over low heating and the beaker reweighed. The concentration of oil and grease is determined by the increase in mass of the collection beaker per volume of sample extracted. O&G is suitable for lubricating oils and other high boiling point products but is not suitable for volatiles. Reference to APHA 5520 B and USEPA 1664 Revision B.. Internal Reference AN185
AN245	Anions by Ion Chromatography: A water sample is injected into an eluent stream that passes through the ion chromatographic system where the anions of interest ie Br, Cl, NO <sub>2</sub> , NO <sub>3</sub> and SO <sub>4</sub> are separated on their relative affinities for the active sites on the column packing material. Changes to the conductivity and the UV-visible absorbance of the eluent enable identification and quantitation of the anions based on their retention time and peak height or area. APHA 4110 B
AN277	Nitrite ions, when reacted with a reagent containing sulphanilamide and N-(1-naphthyl)-ethylenediamine dihydrochloride produce a highly coloured azo dye that is measured photometrically at 540nm.
AN279/AN293(Sydney)	The sample is digested with Sulphuric acid, K <sub>2</sub> SO <sub>4</sub> and CuSO <sub>4</sub> . All forms of phosphorus are converted into orthophosphate. The digest is cooled and placed on the discrete analyser for colorimetric analysis.
AN281	An unfiltered water or soil sample is first digested in a block digester with sulfuric acid, K <sub>2</sub> SO <sub>4</sub> and CuSO <sub>4</sub> . The ammonia produced following digestion is then measured colourimetrically using the Discrete Analyser. A portion of the digested sample is buffered to an alkaline pH, and interfering cations are complexed. The ammonia then reacts with salicylate and hypochlorite to give a blue colour whose absorbance is measured at 660nm and compared with calibration standards. This is proportional to the concentration of Total Kjeldahl Nitrogen in the original sample.

FOOTNOTES

IS	Insufficient sample for analysis.	LOR	Limit of Reporting
LNR	Sample listed, but not received.	↑↓	Raised or Lowered Limit of Reporting
*	NATA accreditation does not cover the performance of this service.	QFH	QC result is above the upper tolerance
**	Indicative data, theoretical holding time exceeded.	QFL	QC result is below the lower tolerance
***	Indicates that both * and ** apply.	-	The sample was not analysed for this analyte
		NVL	Not Validated

Unless it is reported that sampling has been performed by SGS, the samples have been analysed as received.  
Solid samples expressed on a dry weight basis.

Where "Total" analyte groups are reported (for example, Total PAHs, Total OC Pesticides) the total will be calculated as the sum of the individual analytes, with those analytes that are reported as <LOR being assumed to be zero. The summed (Total) limit of reporting is calculated by summing the individual analyte LORs and dividing by two. For example, where 16 individual analytes are being summed and each has an LOR of 0.1 mg/kg, the "Totals" LOR will be 1.6 / 2 (0.8 mg/kg). Where only 2 analytes are being summed, the "Total" LOR will be the sum of those two LORs.

Some totals may not appear to add up because the total is rounded after adding up the raw values.

If reported, measurement uncertainty follow the ± sign after the analytical result and is expressed as the expanded uncertainty calculated using a coverage factor of 2, providing a level of confidence of approximately 95%, unless stated otherwise in the comments section of this report.

Results reported for samples tested under test methods with codes starting with ARS-SOP, radionuclide or gross radioactivity concentrations are expressed in becquerel (Bq) per unit of mass or volume or per wipe as stated on the report. Becquerel is the SI unit for activity and equals one nuclear transformation per second.

Note that in terms of units of radioactivity:

- a. 1 Bq is equivalent to 27 pCi
- b. 37 MBq is equivalent to 1 mCi

For results reported for samples tested under test methods with codes starting with ARS-SOP, less than (<) values indicate the detection limit for each radionuclide or parameter for the measurement system used. The respective detection limits have been calculated in accordance with ISO 11929.

The QC and MU criteria are subject to internal review according to the SGS QAQC plan and may be provided on request or alternatively can be found here: [www.sgs.com.au/en-gb/environment-health-and-safety](http://www.sgs.com.au/en-gb/environment-health-and-safety).

This document is issued by the Company under its General Conditions of Service accessible at [www.sgs.com/en/Terms-and-Conditions.aspx](http://www.sgs.com/en/Terms-and-Conditions.aspx). Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein.

Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client only. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law .

This report must not be reproduced, except in full.